


Gravimetric Analysis

Definition of Gravimetric Analysis :

It is a type of quantitative analysis that involves weighing of the constituent under determination. IN OTHER WORDS, It is the process of isolating and weighing an element or compound in a pure form.



⇒ Remember:

Classification of Analytical Chemistry




Comparison between

Volumetric Analysis (Titration) & Gravimetric Analysis

	Volumetric Analysis (Titration)	Gravimetric Analysis
Diagram:		
Sample:	Liquid or Solid	Liquid or Solid
Result obtained:	E.P. (Vol. consumed from titrant) Ex: 4 ml	Weight of the ppt Ex: 4 g
Required:	Concentration of the sample	Concentration of the sample
Steps of Calculations:	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1st</div> : Calc. of equivalence factor (F). <div style="border: 1px solid black; padding: 2px; display: inline-block;">2nd</div> : Calc. of Concentration.	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1st</div> : Calc. of gravimetric factor (F). <div style="border: 1px solid black; padding: 2px; display: inline-block;">2nd</div> : Calc. of Concentration.

Examples of gravimetric determinations :


1) Gravimetric determination of NaCl:



ppt

↓
filtered, dried and weighed and then from the weight of the ppt, the concentration of NaCl sample can be calculated.

2) Gravimetric determination of Potash alum [KAl(SO₄)₂.12H₂O]:



gelatinous ppt

↓
Not suitable for weighing.

ppt

↓
More suitable for weighing & from its weight, the concentration of Potash alum sample can be calculated.


Calculations for Gravimetric Analysis { 2 STEPS }:

1st step: Calculation of the Gravimetric Factor OR Conversion Factor (F):

Gravimetric Factor is the weight of the sample equivalent 1 g of the ppt.


$$F = \frac{m \times (\text{M.W. of substance sought})}{n \times (\text{M.W. of substance weighed})}$$

sample
final ppt



m & n are integers that make the molecular weights in the numerator (البسط) & denominator (المقام) chemically equivalent.

N.B. Gravimetric factor (F) is completely different from Equivalence factor (F).



2nd step: Calculation of the Concentration:

If the sample is liquid (solution):

$$\text{Concn.} = \frac{\text{Wt. of the ppt} \times F \times 1000}{\text{Volume taken from the sample}} = \dots\dots\dots \text{g/L}$$


$$= \frac{\text{Wt. of the ppt} \times F \times 100}{\text{Volume taken from the sample}} = \dots\dots\dots \text{g\% (i.e. g/100 ml)}$$

If the sample is solid (powder):

$$\text{Concn. or Purity} = \frac{\text{Wt. of the ppt} \times F \times 100}{\text{Wt. taken from the sample}} = \dots\dots\dots \text{g\% (i.e. g/100 g)}$$

Example: In the determination of $MgSO_4 \cdot 7H_2O$:

The reaction can be summarized as follows:



$$\begin{array}{ccc}
 2 \text{ (MgSO}_4 \cdot 7\text{H}_2\text{O)} & \equiv & 1 \text{ (Mg}_2\text{P}_2\text{O}_7) \\
 \downarrow & & \downarrow \\
 \text{m} & & \text{n} \\
 \text{sample} & & \text{final ppt}
 \end{array}$$


$$F = \frac{2 \times (\text{M.W. of MgSO}_4 \cdot 7\text{H}_2\text{O})}{1 \times (\text{M.W. of Mg}_2\text{P}_2\text{O}_7)}$$

Then the concentration can be calculated from one of the above mentioned equations:

$$\text{Concn.} = \frac{\text{Wt. of the ppt} \times F \times 1000}{\text{Volume taken from the sample}} = \dots\dots\dots \text{ g/L}$$




Derivation of Gravimetric Factor (F):

Example: In the determination of $MgSO_4 \cdot 7H_2O$:



The Gravimetric Factor (F) means:

the weight of the sample ($MgSO_4 \cdot 7H_2O$) equivalent 1 g of the ppt ($Mg_2P_2O_7$).



$$\therefore F = 1 \times \frac{2 \times (\text{M.W. of MgSO}_4 \cdot 7\text{H}_2\text{O})}{1 \times (\text{M.W. of Mg}_2\text{P}_2\text{O}_7)}$$

$$\therefore F = \frac{2 \times (\text{M.W. of MgSO}_4 \cdot 7\text{H}_2\text{O})}{1 \times (\text{M.W. of Mg}_2\text{P}_2\text{O}_7)}$$

→ sample
→ final ppt

Basic Steps for Any Gravimetric Analysis

1. Sampling & Preparation for precipitation.
2. Precipitation.
3. Digestion (i.e. Leaving the ppt in contact with the soln. for a considerable time).
4. Filtration & Washing.
5. Drying OR Ignition.
6. Weighing & Calculations.

Required Glassware for Gravimetric Analysis (الأدوات المطلوبة):


- Beaker 250 ml.
- 10-ml bulb pipette.
- 10-ml graduated pipette.
- 10-ml measuring cylinder.
- Filter papers (Whatmann no.1).
- Funnel.
- Flask.
- Watch glass.
- Test tube.

Exp. (8): Gravimetric Determination of $\text{Na}_2\text{SO}_4 \cdot 10 \text{H}_2\text{O}$

(USP 2007)

Principle :

It depends on the precipitation of SO_4^{--} as BaSO_4 by adding BaCl_2 to a hot solution of SO_4^{--} slightly acidified with HCl .



Procedure :

Sampling &
Preparation
for
Precipitation

1- Transfer 10 ml of the sample of sodium sulphate into a 250 ml beaker, Dilute to about 150 ml with water.

2- Add about 10 ml dilute HCl , and heat to boiling.

3- Remove the flame, and add slowly with constant stirring (*by glass rod*), a hot 2 % BaCl_2 solution in slight excess (*about 10 ml heated in a test tube*).

Precipitation

4- **Test for complete precipitation** by leaving the precipitate to settle, and then adding to the clear supernatant liquid, drops of the precipitating agent (2% BaCl_2), noting if any turbidity appears. *If turbidity appears* → add more drops of BaCl_2 soln. and repeat the test.

If no turbidity → complete the next steps.



Digestion

5- Cover the beaker with a watch glass, and place it in a water bath almost at the boiling point for about one hour.

Filtration & Washing

6- Decant the supernatant liquid in a suitable ashless filter paper and wash the precipitate in the beaker by decantation two times with hot water.


Decantation means: *Pouring off the clear upper portion (supernatant) of a fluid leaving a sediment or a ppt at the bottom*
WHILE Washing by Decantation means: *mixing some of the washing solution with the ppt, then allow the ppt to settle down again and then the clear supernatant is decanted.*

Repeat washing & filtration until all the ppt is filtered

Drying

Weighing & Calculations

Calculations :



$$F = \frac{1 \times (\text{M.W. of Na}_2\text{SO}_4 \cdot 10\text{H}_2\text{O})}{1 \times (\text{BaSO}_4)}$$

$$\text{Concn.} = \frac{\text{Wt. of the ppt} \times F \times 1000}{\text{Volume taken from the sample}} = \dots\dots\dots \text{ g/L}$$

***** >>Best wishes >> *****